

orthocare®



orthopaedics

CONTENT

Features 03

Materials 04 Air-X
05 Comfort / Easy Knitted
06 Comfort Polypropylene
07 Comfort Mycrofiber


New 08

Products 18 Premium Line
20 Kids
22 Basic and Sports
35 Pillows




Soft touch

Knitting technology and special materials increases wearing comfort and provide greater feeling.



Air flow

Prevents itches and irritations on the skin by the way of active breathing.




Compression

Makes swellings are relieved and the pain is reduced together with proprioceptive effect.



Perfect fit

Product fits anatomically and provides optimum medical care.



Variety

Different solutions for different problems to reply all needs



Comfort Zone

Increases comfort and satisfaction thanks to specially knitted zones.

AirX

Let your body breathe!

AirX is an innovative, latex-free alternative to neoprene, foam and elastics. AirX is highly breathable, and comfortable to the skin, leading to increased patient comfort and compliance.

AirX offers the support without the issues of moisture build up. The ideal solution for the orthopaedic and sports medicine.

Features:


- 4-way stretch
- High breathable
- Effective moisture vapour transport
- Low weight
- Latex free
- Heat mouldable
- Rigid surface
- Washable



Comfort & Easy Knitted

Breathing knitted technology!


- Breathable
- Moisture transport
- Highly elastic
- Perfect compliance due to specially knitted zones
- Massage effect of compressive fabric
- Easy to wear and use
- Comfortable
- Anatomically knitted.



Comfort Polypropylene

Smart fabric!


- Vapor transport
- High breathable
- Lightweight
- Anti-allergical
- Anti-electrostatic
- Insensitive to absorb liquids (transfers moisture)
- Anti-smell
- Resistance to light
- Comfortable



Comfort Microfiber

Breathing knitted technology!

- Breathable: provide moisture vapour transport
- Skin friendly: much more comfortable than all other stockinettes
- Provide high patient compliance
- Durable and elegant



Genucare Air-X ROM Genucare Air-X ROM wrap

Knee brace

Extension limiting wedges:

0°, 10°, 20°, 30°, 45°

Flexion limiting wedges:

0°, 30°, 45°, 60°, 75°, 90°

Features

- breatheable and moisture transferring innovative Air-X material
- more compliance with comfort zone at popliteal area
- stabilisation provided by medial and lateral ROM hinges and wrap around loop straps.
- 6775: easy to fit and remove thanks to anatomically shaped wrap design
- 6770: sleeve design



NEW

Walker boot short Walker boot long

Lower-leg foot brace

Features

- Easy to use due to adjustable loop straps
- Protection and comfort by padding
- Physiological gait thanks to rocker sole design
- Range of Motion hinge (7455)
- 4 sizes

ROM Walker

Lower-leg foot brace with ROM hinge

Plantarflexion:

0°, 10°, 20°, 30°, 40°

Dorsiflexion:

20°, 10°, 0°

Immobilisation:

at each 10°




Malleocare Pro

Ankle brace

Features

- Lateral and medial stabilizing rods guarding malleolus
- Lacing enabling individual adjustment of the brace in 2 colours black and magenta.
- The two side strap system stabilising the ankle
- Anatomical design and padding for comfort and compliance.
- High quality durable loop fabric and straps.



Dorsocare Inno

Lumbar + Spinal support

Features:

- 2 in 1 new and fancy design
- Innovative production process
- Advanced accessories such as hand grips
- Breathable, premium quality of elastics
- Hook resistant fabrics
- Easy to fit and remove
- Higher stabilisation with adjustable straps
- Supported with flexible ribs and steel stays.
- Detachable, free use as lumbar or spinal support

2 pieces set



NEW

Manucare Easy, Epicare Easy, Malleocare Easy, Genucare Easy

Bright and fancy colored elastic supports.

- Breathable
- Highly elastic
- Seamless
- Perfect compliance due to specially knitted zones
- Massage effect of compressive fabric
- Easy to wear and use.



Manucare III stable

Wrist support with stabilizing strap

Manucare III plus

Wrist support with thumb

Features

- Easy to fit and remove sleeve design
- High compliance and comfort due to anatomical shape
- Practical velcro system and loops for pressure adjustment
- Wrist strap for more stabilisation (4542)



Malleocare step

Ankle brace

Features

- Easy to use due to adjustable velcro straps
- Perfect wearing comfort due to padding.
- Anatomically shaped protective shells at lateral and medial.
- Universal design that fits to all sizes.



Dosicare maternity

Dosicare maternity plus

Lumbar supports for pregnancy

Features

- Cushioning effect with high quality soft fabric
- Easy to fit and remove due to practical velcro system and hand loops
- Higher stabilisation with adjustable upper straps
- Flexible ribs for slight limitation.
- Adjustable groin straps for more comfort (Model.2720)




PREMIUM





1170 Reinforced collar

Reinforced cervical collar



Indications

Cervical syndrome (degenerative, posttraumatic) cervical disc hernia, cervical strain/spasm, tension headache.

Mode of Action

Supports cervical joints and muscles, provides heating and proprioceptive effect, slightly restricts cervical movements in all directions.


Features

Air permeable microfiber fabric, high density sponge, polyethylene splint. In colors beige and black.



1180 Soft collar

Soft sponge cervical collar



Indications

Cervical syndrome (degenerative, posttraumatic) cervical disc hernia, cervical strain/spasm, tension headache.

Mode of Action

Supports soft tissues of cervical region, provides heating and proprioceptive effect, slightly restricts cervical movements in all directions.


Features

Air permeable microfiber fabric, high density sponge. In colors beige and black.



1330 Clavicare

Clavicula support



Indications

Clavicular fractures, on shoulder joint injuries.

Mode of Action

Keeps shoulders in normal alignment and helps for healing clavicular fracture, supports thoracic spine in normal position.


Features

Comfortable and adjustable straps, high quality and durable terry cloth material, easy to fit and remove due to practical velcro system, comfortable when lied on back.



1741 Velpocare

Shoulder immobiliser



Indications

Humerus fractures, shoulder stabilisation.

Mode of Action


Stabilises shoulder complex, prevents inferior subluxation/dislocation of glenohumeral joint.

Features

Comfortable and adjustable straps, high quality and durable terry cloth material, easy to fit and remove due to practical velcro system, additional stabilisation by the chest strap.



3510 Armsling velcro



Indications

Sub-acromial decompressions, global shoulder instability, rotator cuff strains, shoulder dislocation/subluxation, shoulder separation, non-displaced humeral fractures.

Mode of Action


Stabilizes shoulder complex, prevents inferior subluxation, dislocation of glenohumeral joint.

Features

Comfortable and adjustable straps, high quality and durable terry cloth material, easy to fit and remove due to practical velcro system, more comfortable with thumb sleeve.



3511 Shoulder abduction 15



Indications

Post-operative immobilisation, post-op rotator cuff reconstruction, subacromial decompression, conservative immobilisation of the shoulder.

Mode of Action

Provides special positioning for the shoulder complex in 15 degrees abduction, prevent subluxation of shoulder.

Features

Comfortable and adjustable straps, high quality and durable terry cloth material, easy to fit and remove due to practical velcro system, portable arm sling, perfect positioning due to shaped cushion.



3512 Shoulder abduction 30/45

universal

Indications

Post-operative immobilisation, post-op rotator cuff reconstruction, subacromial decompression, conservative immobilisation of the shoulder.

Mode of Action

Provides special positioning for the shoulder complex in 30 or 45 degrees abduction, prevent subluxation of shoulder.

Features

Comfortable and adjustable straps, high quality and durable terry cloth material, easy to fit and remove due to practical velcro system, portable arm sling, optional positioning at 30 or 45 degrees.



3513 Shoulder abduction lateral

universal

Indications

Post-operative immobilisation, post-op rotator cuff reconstruction, subacromial decompression, conservative immobilisation of the shoulder.

Mode of Action

Provides special positioning for the shoulder complex in 90 degrees abduction at laterals,

prevents subluxation of shoulder.

Features

Comfortable and adjustable straps, high quality and durable terry cloth material, easy to fit and remove due to practical velcro system, portable arm sling, perfect positioning due to shaped cushion.



1335 Acromio clavicare

Acromioclavicular brace

universal

Indications

Clavicle fractures, on shoulder joint injuries.

Mode of Action

Keeps shoulders in normal alignment and helps for healing clavicular fracture, supports thoracic spine in normal position.

Features

Comfortable and adjustable straps, high quality and durable terry cloth material, easy to fit and remove due to practical velcro system.



3910 Epicare Comfort

Knitted elbow support with silicon pads

cm	
S	: 20-23
M	: 23-27
L	: 27-30
XL	: 30-32
XXL	: 32-35

Indications

Tenosynovitis, osteoarthritis/arthritis, elbow pain, medial and lateral epicondylitis

Mode of Action

Helps for maintaining correct joint position and muscle support during the function, supports soft tissues around

the elbow, provides proprioceptive effect and healing process.

Features

Anatomically knitted and supported with lateral pads, specially knitted zones at the bend of elbow, high elasticity, breathable comfort knitting system.



3920 Epicare Comfort plus

Knitted elbow support with silicon pads and epicondylitis strap

cm	
S	: 20-23
M	: 23-27
L	: 27-30
XL	: 30-32
XXL	: 32-35

proprioceptive effect and healing process.

Features

Anatomically knitted and supported with lateral pads, specially knitted zones at the bend of elbow, removable retaining strap for pressure adjustment by velcro system, high elasticity, breathable comfort knitting system.



3110 Epicare brace

Elbow brace

universal

Indications

Medial and lateral epicondylitis, epicondylitis humeri (e.g. tennis elbow).

Mode of Action

Helps for maintaining correct joint position and muscle support during the function, supports soft tissues around the elbow, provides proprioceptive effect and

healing process by compression.

Features

Innovative, simple and effective design, easy to fit, adjustable pressure by loop system and portable pad.



3120 Epicare

Elbow strap with single pad

universal

Indications

Medial and lateral epicondylitis, epicondylitis humeri (e.g. tennis elbow).

Mode of Action

Helps for maintaining correct joint position and muscle support during the function, supports soft tissues around the elbow, provides proprioceptive effect.


Features

Simple, easy to fit and effective design, adjustable pressure by velcro system and removable pad.



4910 Manucare II Comfort

Knitted wrist support with stabilizing strap

cm	left/right
	S : 14-16 M : 16-18 L : 18-20 XL : 20-22

Indications

Tenosynovitis, carpal instabilities, CTS, after casting, arthralgia (reactive, rheumatic) ganglion pain, traumas of soft tissues of wrist.

Mode of Action

Provides more stabilisation with extra strap, free range of motion for the thumb.


Features

Breathable comfort knitting system, high elasticity, anatomically shaped and fitted design, high compliance with specially knitted zones, wide strap for pressure adjustment by velcro system reformable thermoplastic splint.



4515 Manucare Comfort

Wrist support

cm	left/right
	XS : 14-15 S : 15-16,5 M : 16,5-18 L : 18-19,5 XL : 19,5-21 XXL : 21-22,5

Indications

Tenosynovitis, carpal instabilities, CTS, after casting, arthralgia (reactive, rheumatic) ganglion pain, traumas of soft tissues of wrist.

Mode of Action

Stabilises the wrist, free range of motion for the thumb.


Features

High compliance and comfort thanks to breathable, transferring moisture and anti-allergical high-tech polypropylene fabric, anatomically shaped and fitted design, reformable thermoplastic splint.



4525 Manucare Comfort plus

Wrist support with thumb

cm	left/right
	XS : 14-15 S : 15-16,5 M : 16,5-18 L : 18-19,5 XL : 19,5-21 XXL : 21-22,5

Indications

MCP joint disorders, CMC joint degenerations, de quarvein, CTS + de quarvein.

Mode of Action

Stabilises the wrist and thumb.


Features

High compliance and comfort thanks to breathable, transferring moisture and anti-allergical high-tech polypropylene fabric, anatomically shaped and fitted design, reformable thermoplastic splint.



4535 Thenarcare Comfort

Thenar support

cm	left/right
	XS : 13-14 S : 14-16 M : 16-18 L : 18-20 XL : 20-22

Indications

Degenerative changes of 1 st MCP, de quarvein, osteoarthritis, traumatic disorders.

Mode of Action

Stabilises the MCP and CMP joints of thumb.


Features

High compliance and comfort thanks to breathable, transferring moisture and anti-allergical high-tech polypropylene fabric, anatomically shaped and fitted design.



4545 Manucare Comfort stable

Wrist support with stabilizing strap and wide thermoplastic splint

cm	left/right
	XS : 14-15 S : 15-16,5 M : 16,5-18 L : 18-19,5 XL : 19,5-21 XXL : 21-22,5

Indications

Tenosynovitis, carpal instabilities, CTS, after casting, arthralgia (reactive, rheumatic) ganglion pain, traumas of soft tissues of wrist.

Mode of Action

Provides more stabilisation with extra strap, free range of motion for the thumb.


Features

High compliance and comfort thanks to breathable, transferring moisture and anti-allergical high-tech polypropylene fabric, anatomically shaped and fitted design, wide strap for pressure adjustment by velcro system, wide and reformable thermoplastic splint.



1320 Posturecare

Posture support

cm	
	S : 60-75 M : 75-90 L : 90-105 XL : 105-120

Indications

As postural support, poor posture, weak thoracic muscles.

Mode of Action

Helps for maintaining up-right posture, provides proprioceptive feed back for back muscles.

Features

Skin friendly, high quality elastic fabric, easy to fit and remove due to practical velcro system, adjustable shoulder straps, soft spring stays for slight limitation.




1710 Sternicare

Sternum support with flexible ribs

1711 Sternicare plus

Sternum support with rods

cm	
	S : 65-75 M : 75-85 L : 85-95 XL : 95-105 XXL : 105-115

Indications


After surgery, when sternum stabilization needed.

Mode of Action

Stabilises sternoclavicular joints and ribcage, helps for healing in post-surgery period by providing heating and proprioceptive effect.


Features

Skin friendly, high quality elastic fabric, easy to fit and remove due to practical velcro system, adjustable shoulder straps.



2550 Dorsocare

Dorsolomber back support with steel splints

cm	
	S : 60-75 M : 75-90 L : 90-105 XL : 105-120 XXL : 120-135

Indications

Thoracic and lumbar stabilisation, postoperative immobilisation, acute and chronic lumbago.

Mode of Action

Supports and keeps spine correct postural position, restricts slightly spinal

movements in all directions and gives proprioceptive feed back.


Features

Skin friendly, high quality elastic fabric, easy to fit and remove due to practical velcro system and hand loops, higher stabilisation with adjustable straps, steel splints for limitation.



2510 Abdominocare

Abdominal support

cm	
	S : 60-75 M : 75-90 L : 90-105 XL : 105-120 XXL : 120-135

Indications

Post-operative.

Mode of Action

Keeps abdominal muscles in inner range and prevents further weakness, support abdominal wall for healing after a surgery.


Features

Skin friendly, high quality elastic fabric, easy to fit and remove due to practical velcro system.



2540 Lumbocare high

High lumbosacral back support with steel splints

cm	
	S : 60-75 M : 75-90 L : 90-105 XL : 105-120 XXL : 120-135

Indication

Accute back pain, degenerative spinal changes, paravertebral spasm, iliosacral irritations, lumbago with insufficiency, tendomyopathies.

Mode of Action

Provides support from pelvis to mid thoracic region, gives feed back for postural alignment and movement patterns.


Features

Skin friendly, high quality elastic fabric, easy to fit and remove due to practical velcro system and hand loops, higher stabilisation with adjustable straps, steel splints for limitation.



2530 Lumbocare light

Lumbosacral back support with soft splints

cm	
	S : 60-75 M : 75-90 L : 90-105 XL : 105-120 XXL : 120-135

Indications

Accute and chronic lumbago, lumbago siciatica, after lumbar traction, hypermobile and degenerative changes.

Mode of Action

Provides support from pelvis to mid thoracic region, gives feed back for postural alignment and movement patterns.


Features

Skin friendly, high quality elastic fabric, easy to fit and remove due to practical velcro system and hand loops, higher stabilisation with adjustable straps, soft springs for slight limitation.



2526 Lumbocare Inno-II

High lumbosacral back support with double tightening

cm	
	S : 60-75 M : 75-90 L : 90-105 XL : 105-120 XXL : 120-135

Indications

Accute back pain, degenerative spinal changes, paravertebral spasm, iliosacral irritations, lumbago with insufficiency, tendomyopathies.

Mode of Action

Provides support from pelvis to mid thoracal region, gives feed back for postural alignment and movement patterns.


Features

Skin friendly, breathing and high quality elasticity thanks to the fine fabric, easy to fit and remove due to practical velcro system and hand loops, higher stabilisation with adjustable straps, tightening thanks to the additional strap, supported with steel stays.



2521 Lumbocare Inno-I

Lumbosacral back support with double tightening

cm	
	S : 60-75 M : 75-90 L : 90-105 XL : 105-120 XXL : 120-135

Indications

Accute and chronic lumbago, lumbago siciatica, after lumbar traction, hypermobile and degenerative changes.

Mode of Action

Provides support from pelvis to mid thoracal region, gives feed back for postural alignment and movement patterns.


Features

Skin friendly, breathing and high quality elasticity thanks to the fine fabric, easy to fit and remove due to practical velcro system and hand loops, higher stabilisation with adjustable straps, tightening thanks to the additional strap, supported with steel stays.



2560 Dorsocare Inno

Lumbar + Spinal support with rods and flexible ribs

cm	
	S : 60-75 M : 75-90 L : 90-105 XL : 105-120 XXL : 120-135



Indications

Thoracic and lumbar stabilisation, postoperative immobilisation, acute and chronic lumbago.

Mode of Action

Supports and keeps spine correct postural position, restricts slightly spinal movements in all directions and gives proprioceptive feed back.


Features: 2 supports in 1 article.

Skin friendly, breathing and high quality elasticity thanks to the fine fabric, easy to fit and remove due to practical velcro system and hand loops, higher stabilisation with adjustable straps, tightening thanks to the additional strap, supported with rods.



2910 Lumbocare Comfort

Lumbosacral back support with splints

cm	
	S : 60-75 M : 75-90 L : 90-105 XL : 105-120 XXL : 120-135

Indications

Acute back pain, degenerative spinal changes, paravertebral spasm, iliosacral irritations, lumbago, postdistectomy syndrome.

Mode of Action

Support lumbar spine and back mussels, provide

proprioceptive effect, slightly reduces lordosis.

Features

Breathable comfort knitting system, high compliance with anatomically knitted design, high elasticity, supported with steel stays, easy to fit and remove due to practical velcro system.



2110 Dasicare

Lumbosacral back support with flexible ribs

2130 Dasicare plus

Lumbosacral back support with rods

cm	
	S : 64-72 M : 72-80 L : 80-90 XL : 90-100 XXL : 100-110

Indications

Acute back pain, degenerative disc, paravertebral muscle spazm, weak back muscles to increase proprioception.

Mode of Action

Supports joints and soft tissues of lumbar region, provides heating and proprioceptive effect, restricts especially

flexion and extension movements in lumbar spine.

Features

Skin friendly, high quality and durable elastic fabric, easy to fit and remove due to practical velcro system, high compliance at abdominal region with divided straps, steel splints for strong limitation.




2710 Dasicare maternity

Lumbar support for pregnancy

2720 Dasicare maternity plus

Lumbar support for pregnancy with groin bands

cm	
	S : 70-85 M : 85-100 L : 100-115

Indications

Acute back pain, degenerative spinal changes, paravertebral spasm, iliosacral irritations, lumbago with insufficiency, tendomyopathies.

Mode of Action

Provides support from pelvis to mid thoracal region, gives feed back for postural alignment and movement patterns.

Features

Cushioning effect with high quality soft fabric, easy to fit and remove due to practical velcro system and hand loops, higher stabilisation with adjustable upper straps, flexible ribs for slight limitation.



6910 Genucare Comfort

Knitted elastic knee support with patella silicon and soft splints

	A	B
XS	41-43	28-31
S	43-46	31-34
M	46-49	34-37
L	49-52	37-40
XL	52-56	40-43
XXL	56-58	43-46

Indications

Arthritis, osteoarthritis, irritation posttraumatic and post-operative, tendomyopathies, feeling of instability.

Mode of Action

Supports and stabilises the kneecap and joint, patella ring and patella cut out give extra

support, comfort material improves proprioception.

Features

Breathable comfort knitting system, high compliance with specially knitted zones at the popliteal area, high elasticity, supported with medial and lateral spiral stays.



6920 Genucare Comfort plus

Knitted elastic knee support with patella silicon, soft splints and stabilizing straps

	A	B
XS	41-43	28-31
S	43-46	31-34
M	46-49	34-37
L	49-52	37-40
XL	52-56	40-43
XXL	56-58	43-46

Indications

Arthritis, osteoarthritis, irritation posttraumatic and post-operative, tendomyopathies, feeling of instability.

Mode of Action

Supports and stabilises the kneecap and joint, patella ring and patella cut out give extra support, straps provide

stabilisation and wearing comfort, improves proprioception.

Features

Breathable comfort knitting system, high compliance with specially knitted zones at the popliteal area, high elasticity, supported with medial and lateral spiral stays.



6730 Genucare AirX ligament

Knee support with patella pad and soft splints

cm	
S	30-34
M	34-38
L	38-42
XL	42-46
XXL	46-50

Indications

Arthritis, osteoarthritis, irritation posttraumatic and post-operative, tendomyopathies, feeling of instability.

Mode of Action

Supports and stabilises the kneecap and joint, patella ring and patella cut out give extra support, Air-X material improves proprioception.

Features

Breathable and moisture transferring innovative elastic Air-X material, high compliance due to wrap around velcro straps, easy to pull on and comfortable wearing due to special popliteal zone material, medial and lateral spiral stays.



6735 Genucare AirX ligament wrap

Knee support with patella pad and soft splints

cm	
S	30-34
M	34-38
L	38-42
XL	42-46
XXL	46-50

Indications

Arthritis, osteoarthritis, irritation posttraumatic and post-operative, tendomyopathies, feeling of instability.

Mode of Action

Supports and stabilises the kneecap and joint, patella ring and patella cut out give extra support, Air-X material improves proprioception.

Features

Breathable and moisture transferring innovative elastic Air-X material, high compliance due to wrap around velcro straps, easy to pull on and comfortable wearing due to special popliteal zone material, medial and lateral spiral stays.



6750 Genucare AirX stable

Knee support with patella pad and hinge

cm	
S	30-33
M	34-38
L	38-42
XL	42-46
XXL	46-50

Indications

ACL and PCL injuries, 1st and 2nd grade meniscus degenerations, collateral ligament injury, arthritis (progressive rheumatoid arthritis).

Mode of Action

Supports and stabilizes the knee and joint, patella ring and patella cut out give extra support, Air-X material

improves proprioception, prevents hyperextension.

Features

Breathable and moisture transferring innovative elastic Air-X material, high compliance due to wrap around velcro straps, easy to pull on and comfortable wearing due to special popliteal zone material, medial and lateral polycentric hinges.



6755 Genucare AirX wrap

Knee support with hinge

cm	
S	30-34
M	34-38
L	38-42
XL	42-46
XXL	46-50

Indications

ACL and PCL injuries, 1st and 2nd grade meniscus degenerations, collateral ligament injury, arthritis (progressive rheumatoid arthritis).

Mode of Action

Supports and stabilises the knee and joint, Air-X material improves proprioception, prevents hyperextension.

Features

Breathable and moisture transferring innovative elastic Air-X material, high compliance due to wrap around version and practical velcro closing system at the front, comfortable wearing due to special popliteal zone material, medial and lateral polycentric hinges.



6770 Genucare AirX ROM

Knee support with patella pad and ROM hinge

- **Extension:** 0°, 10°, 20°, 30°, 45°
- **Flexion:** 0°, 30°, 45°, 60°, 75°, 90°

cm	S : 30-34
	M : 34-38
	L : 38-42
	XL : 42-46
	XXL : 46-50

Indications

ACL and PCL injuries, 1st and 2nd grade meniscus degeneration, collateral ligament injury, arthritis (progressive rheumatoid arthritis), meniscus injuries, in traumatic or degenerative knee joint changes, mild-moderate osteoarthritis

Mode of Action

Provides secure fit at the knee and joint with of extension and flexion limitation, support and stabilizes knee, collateral

stabilisation with sidebars, patella ring and patella cut out give extra support, Air-X material improves proprioception, prevents hyperextension.

Features

Breatheable and moisture transferring innovative Air-X material, easy to fit and remove, thanks to anatomically shaped wrap design, more compliance with comfort zone at popliteal area and wrap around loop straps, stabilisation provided by comprehensive set of wedges for extension and flexion limitation.



6775 Genucare AirX ROM wrap

Knee support with ROM hinge

- **Extension:** 0°, 10°, 20°, 30°, 45°
- **Flexion:** 0°, 30°, 45°, 60°, 75°, 90°

cm	S : 30-34
	M : 34-38
	L : 38-42
	XL : 42-46
	XXL : 46-50

Indications

ACL and PCL injuries, 1st and 2nd grade meniscus degeneration, collateral ligament injury, arthritis (progressive rheumatoid arthritis), meniscus injuries, in traumatic or degenerative knee joint changes, mild-moderate osteoarthritis

Mode of Action

Provides secure fit at the knee and joint with of extension and flexion limitation, support and stabilizes knee, collateral

stabilisation with sidebars, Air-X proprioception, prevents hyperextension.

Features

Breatheable and moisture transferring innovative Air-X material, easy to fit and remove thanks to anatomically shaped wrap design, more compliance with comfort zone at popliteal area and wrap around loop straps, stabilisation provided by comprehensive set of wedges for extension and flexion limitation.



6741 Genucare AirX patella

Knee support with patella pad

cm	S : 30-34
	M : 34-38
	L : 38-42
	XL : 42-46
	XXL : 46-50

Indications

Osteoarthritis, effusion, other painful conditions of knee.

Mode of Action

Knee support with pad that gives extra support and stability to the kneecap and joint, warmth, improves proprioception.

Features

Breatheable and moisture transferring innovative elastic Air-X material, anatomically shaped and fitted.



6760 Genucare AirX hyperX

Knee brace for prevention of hyper extension

cm	S : 30-33
	M : 33-36
	L : 36-39
	XL : 39-42
	XXL : 42-45

Indications

ACL and PCL injuries, 1st and 2nd grade meniscus degeneration, collateral ligament injury, arthritis (progressive rheumatoid arthritis).

Mode of Action

Supports and stabilises the knee and joint, patella ring and patella cut out give extra support, Air-X material improves proprioception,

prevents hyperextension.

Features

Breatheable and moisture transferring innovative elastic Air-X material, prevention of hyperextension thanks to crossed straps at the back, high compliance due to velcro straps at the front, easy to pull on and comfortable wearing due to special popliteal zone material, medial and lateral polycentric hinges.



6310 Genucare immobiliser

Knee immobiliser

cm	S : 40
	M : 50
	L : 60

Indications

Postoperative immobilisation, meniscus injuries, ACL and PCL management.

Mode of Action

Immobilizes the knee, provides less pressure on the patello femoral joint.

Features

Easy adjustment of wrap around straps due to practical velcro system, removable medial, lateral and posterior stabilisation rods, high compliance due to skin-friendly and durable terry cloth material.



6325 Genucare ROM-T

Knee brace with flexion/extension limitations and telescopic length adjustment

Extension Limitation:
-15°, 0°, 15°, 30°, 45°, 60°, 75°
Flexion Limitation:
0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°
universal

Indications

Post operative immobilisation, immobilisation after knee injuries, secure stabilisation of injured ligaments.

Mode of Action

Immobilisation, controlled mobilisation, collateral stabilisation of knee joint.

Features

Skin friendly padding, easy to adjust hinge, easy to adjust length thanks to telescopic system, wide range of flexion and extension limitations.



7910 Malleocare Comfort

Knitted ankle support with silicon pads

cm	
XS	: 17-20
S	: 20-22
M	: 22-24
L	: 24-26
XL	: 26-28
XXL	: 28-30

Features

Breathable clima-comfort knitting system, anatomically knitted, high compliance with specially knitted zones at the ankle bending area, supported with medial and lateral silicon pads.

Indications

Chronic instability of ankle, ligament injuries, joint effusion/swelling, malleolar sensitivity.

Mode of Action

The pad have a massage effect which helps circulation and increases the blood flow, reduces stiffness and stabilise the ankle, comfort material improves proprioception.



7920 Achicare Comfort

Knitted achilles tendon support with silicon pads

cm	
XS	: 17-20
S	: 20-22
M	: 22-24
L	: 24-26
XL	: 26-28
XXL	: 28-30

Indications

Chronic and acute inflammations of achilles tendon, posttraumatic and post-operative inflammation, achillodynia, achillobursitis.

Mode of Action

The pad have a massage effect which helps circulation and increases the blood flow, reduces stiffness, protects

achilles tendons and stabilises the ankle, comfort material improves proprioception.

Features

Breathable clima-comfort knitting system, anatomically knitted, high compliance with specially knitted zones at the ankle bending area, supported with silicon pad.



7930 Malleocare strap

Knitted ankle support with silicon pads and stabilizing strap

cm	
XS	: 17-20
S	: 20-22
M	: 22-24
L	: 24-26
XL	: 26-28
XXL	: 28-30

more stability to ankle, comfort material improves proprioception.

Features

Breathable clima-comfort knitting, elastic and light material, and anatomically-shaped knit, anatomically-shaped silicone pads for effective pressure relief and stabilisation at the ankle bending area, added stability thanks to wide strap, increased comfort by using softer material in the flexion area.

Indications

Ankle sprains, instability of the ankle joint, joint effusions and swellings, arthritis and osteoarthritis, chronic, post-operative and post-traumatic inflammation

Mode of Action

The pads have a massage effect which helps circulation and increases the blood flow, reduces stiffness and wide strap gives



7210 Foot lifter

Dorsiflexion ankle support

universal

Indications

Foot drop, peroneal neuropathy, needs of dorsiflexion

Mode of Action

Raise the foot to the right position, supports and helps to move.

Features

Simple and effective design, comfortable and easy to use due to elastic velcro straps, high quality material.



7510 Malleocare pro

Ankle brace

cm	
XS	: 18-21
S	: 21-24
M	: 24-27
L	: 27-30

Features

Lateral and medial stabilizing rods guarding malleolus, lacing enabling individual adjustment of the brace, the two-side strap system stabilising the ankle, anatomical design and padding for comfort and compliance

Indications

Chronic ankle instability, ankle sprains, protection of joint after trauma, inflammation of the ankle joint, mild ligament injuries

Mode of Action

Provides very good stabilisation in pronation and supination movements, protects the ankle



7190 Hallux valgus light

Hallux valgus day strap

7191 Hallux valgus combo

Hallux valgus night strap

left/right
universal

Indications

Hallux valgus - conservative and post-operative treatment.

Mode of Action

Strap unload the joints of big toe with a tension and helps moving the big toe to the right position.

Features

Simple and effective design, comfortable and easy to use due to elastic velcro straps, light as to be used in shoes.

7191: Longer straps for more stabilisation during the night time.



7440 Walker boot short

Ankle foot orthosis

shoe size	
S	: 36-39
M	: 40-44
L	: 44-46,5

Indications

Post operative and conservative care after ligament injuries, chronic instabilities of ankle, after distorsions

Mode of Action

Provides stabilisation in pronation and supination movements of ankle, perfect wearing comfort due to padding.

Features

Easy to use due to adjustable velcro straps, protection and comfort by padding.



7445 Walker boot long

Ankle foot orthosis

shoe size	
S	: 36-39
M	: 40-44
L	: 44-46,5

Indications

Chronic and acute inflammations of achilles tendon, posttraumatic and post-operative inflammation, achillobursitis.

Mode of Action

The pad have a massage effect which helps circulation and increases the blood flow, reduces stiffness, protects

achilles tendons and stabilises the ankle, comfort material improves proprioception.

Features

Breathable clima-comfort knitting system, anatomically knitted, high compliance with specially knitted zones at the ankle bending area, supported with silicon pad.



7446 ROM Walker

Ankle foot orthosis

shoe size	
S	: 36-39
M	: 40-44
L	: 44-46,5

Indications

Chronic and acute inflammations of achilles tendon, posttraumatic and post-operative inflammation, achillobursitis.

Mode of Action

The pad have a massage effect which helps circulation and increases the blood flow, reduces stiffness, protects

achilles tendons and stabilises the ankle, comfort material improves proprioception.

Features

Breathable clima-comfort knitting system, anatomically knitted, high compliance with specially knitted zones at the ankle bending area, supported with silicon pad.



7420 Malleocare step

Ankle brace

universal

Indications

Post operative and conservative care after ligament injuries, chronic instabilities of ankle, after distorsions

Mode of Action

Provides stabilisation in pronation and supination movements of ankle, perfect wearing comfort due to padding.

Features

Easy to use due to adjustable velcro straps, protection and comfort by padding.



7410 Foot guard

Ankle brace

shoe size	
S	: 34-37
M	: 38-41
L	: 41-45

Indications

Conditions requiring foot positioning or at risk for toe ulcers, chronic instabilities of ankle

Mode of Action

Easy to adjust T stabilizer helps hold the foot in desired position, extended foot splint helps keep bedding off the toes.

Features

Easy to use due to adjustable velcro straps, protection and comfort by padding, easy to clean removable padding system.



7430 AFO

Ankle foot orthosis

shoe size	
	left/right universal

comfortable due to semi rigid thin design

Features

Lightweight anatomically preshaped, easy to clean, fits in shoes.



KIDS





1339 Clavicare kids

Clavicula support

age	S	: 3-6
	M	: 6-10

Indications

Clavicular fractures, on shoulder joint injuries.

Mode of Action

Keeps shoulders in normal alignment and helps for healing clavicular fracture, supports thoracic spine in normal position.

Features

Comfortable and adjustable straps, high quality and durable terry cloth material, easy to fit and remove due to practical velcro system, comfortable when lied on back.



1749 Velpocare kids

Shoulder immobiliser

age	S	: 3-6
	M	: 6-10

Indications

Humerus fractures, shoulder stabilisation.

Mode of Action

Stabilises shoulder complex, prevents inferior subluxation/dislocation of glenohumeral joint.

Features

Comfortable and adjustable straps, high quality and durable terry cloth material, easy to fit and remove due to practical velcro system, additional stabilisation by the chest strap.



3519 Armsling kids

Arm sling

age	S	: 3-6
	M	: 6-10

Indications

Sub-acromial decompressions, global shoulder instability, rotator cuff strains, shoulder dislocation/subluxation, shoulder separation, non-displaced humeral fractures.

Mode of Action

Stabilizes shoulder complex, prevents inferior subluxation, dislocation of glenohumeral joint.


Features

Comfortable and adjustable straps, high quality and durable terry cloth material, easy to fit and remove due to practical velcro system, more comfortable with thumb sleeve.



4519 Manucare kids

Knitted wrist support with stabilizing strap

cm	left/right
	S : 10-11,5
	M : 11,5-13
	L : 13-14,5

Indications

Tenosynovitis, carpal instabilities, CTS, after casting, arthralgia (reactive, rheumatic) gonglion pain, traumas of soft tissues of wrist.

Mode of Action

Provides more stabilisation with extra strap, free range of motion for the thumb.


Features

Breathable comfort knitting system, high elasticity, anatomically shaped and fitted design, high compliance with specially knitted zones, wide strap for pressure adjustment by velcro system reformable thermoplastic splint.



6319 Genucare immobiliser kids

Knee immobiliser

cm	
	S : 20
	M : 26
	L : 32

Indications

Postoperative immobilisation, meniscus injuries, ACL and PCL management.

Mode of Action

Immobilizes the knee, provides less pressure on the patello femoral joint.

Features

Easy adjustment of wrap around straps due to practical velcro system, removable medial, lateral and posterior stabilisation rods, high compliance due to skin-friendly and durable terry cloth material.



5010 Pavlic strap Ida

month	S	: 2-6
	M	: 4-8
	L	: 6-9

Indications

Harness, congenital hip subluxation and/or dislocation, hip abduction dislocation.

Mode of Action

Keeps hip in flex and abducted position in order to centralize femur against acetabulum.

Features

Comfortable and adjustable straps, high quality and durable material, easy to fit and remove due to practical velcro system, more comfortable with foot place.

BASIC SPORTS





1190 Nelson collar

Nelson cervical collar

		A	B
	0 :	8-10	30-40
	1 :	8-10	35-45
	2 :	8-10	40-50
	3 :	8-10	45-55

Indications

Cervical disc hernia, cervical strain, post-operative, structural loosening of the servical spin.

Mode of Action

Supports cervical joints and muscles, provides heating and proprioceptive effect, slightly restricts cervical movements especially felxion and extension.

Features

X-ray luscent plastozote foam, anatomically shaped



1160 Philadelphia collar

Philadelphia cervical collar

	0 :	4
	1 :	6
	2 :	8
	3 :	10

Indications

Cervical traumas, pre and post-operative, for cervical immobilisation and stabilisation, cervical soft tissue injuries.

Mode of Action

Stabilizes cervical spine, provide heating effect, support soft tissues of cervical region.

Features

X-ray luscent plastozote foam, anatomically shaped back and front parts.



1150 Vitrafix collar

Cervical collar with chin support

		A	B
	1 :	6-8	30-40
	2 :	6-8	35-45
	3 :	8-10	40-50
	4 :	8-10	45-55

Indications

Cervical strains/spazms, cervical disk hernia, torticollis.

Mode of Action

Supports cervical joints and muscles, slightly restricts cervical movements especially at flexion and extension.

Features

High density polyethylene, surrounding with cushioning and skin friendly leather.



1140 Vitra collar

Cervical collar

		A	B
	1 :	6-8	30-40
	2 :	6-8	35-45
	3 :	8-10	40-50
	4 :	8-10	45-55

Indications

Cervical strains/spazms, cervical disk hernia, torticollis.

Mode of Action

Supports cervical joints and muscles, slightly restricts cervical movements especially at flexion and extension.

Features

High density polyethylene, surrounding with cushioning and skin friendly leather.



3410 Armsling

cm		
	S :	34
	M :	38
	L :	44

Indications

Humerus fractures, shoulder injuries, stabilisation of arm, in case of casting.

Mode of Action

Stabilises shoulder complex, prevents inferior subluxation/dislocation of glenohumeral joint.

Features

Durable and adjustable straps, lightweight fabric.



3411 Armsling light

cm		
	S :	30-36
	M :	36-42
	L :	42-48

Indications

Humerus fractures, shoulder injuries, stabilisation of arm, in case of casting.

Mode of Action

Stabilises shoulder complex, prevents inferior subluxation/dislocation of glenohumeral joint.


Features

Durable and adjustable straps, lightweight breathable and comfort fabric.



2610 Dosaicare active

Thermal back support

cm	
	S : 64-72 M : 72-80 L : 80-90 XL : 90-100 XXL : 100-110

Indications

Decreased body stabilisation, insufficient blood circulation, to keep warm.

Mode of Action

Restricts spinal movements in lumbar region, encourages active life style in acute phase, provides support and warmth, increases the blood flow and improves proprioception.


Features

Anatomically shaped and fitted design, high quality neoprene material.



2810 Back support wool

Lumbar support

cm	
	S : 60-75 M : 75-90 L : 90-105 XL : 105-120 XXL : 120-135

Indications

Rheumatism, sciatica, post-traumatic

Mode of Action

Increases the blood flow and eases pain, provides warmth, relieves proprioception.


Features

High quality and elastic angoret.



3130 Epicare active

Elbow support

cm	
	S : 20-23 M : 23-27 L : 27-30 XL : 30-32 XXL : 32-35

Indications

Tenosynovitis, osteoarthritis/arthritis, elbow pain, medial and lateral epicondylitis.

Mode of Action

Helps for maintaining correct joint position and muscle support during the function,

supports soft tissues around the elbow, provides proprioceptive effect and healing process by warming.


Features

Anatomically shaped, simple and comfortable design, high quality of neoprene material.



3135 Epicare active plus

Elbow support with epicondylitis strap

cm	
	S : 20-23 M : 23-27 L : 27-30 XL : 30-32 XXL : 32-35

Indications

Tenosynovitis, osteoarthritis/arthritis, elbow pain, medial and lateral epicondylitis, epicondylitis humeri (e.g. tennis elbow).

Mode of Action

Helps for maintaining correct joint position and muscle support during the function, supports soft tissues around

the elbow, provides proprioceptive effect and healing process by warming.


Features

Anatomically shaped, removable retaining strap for pressure adjustment by velcro system, removable pads, adjustable upper arm width for better compliance.



3911 Epicare Easy

Knitted elastic elbow support

cm	
	S : 20-23 M : 23-26 L : 26-29 XL : 29-32 XXL : 32-35

Indications

Tenosynovitis, osteoarthritis/arthritis, elbow pain, medial and lateral epicondylitis.

Mode of Action

Helps for maintaining correct joint position and muscle support during the function, supports soft tissues around

the elbow, provides proprioceptive effect and healing process.

Features

Simple and comfortable design, specially knitted zones at the bend of elbow, high elasticity, breathable knitting system.



3125 Epicare plus

Elbow strap with double pad

universal

Indications

Medial and lateral epicondylitis, epicondylitis humeri (e.g. tennis elbow).

Mode of Action

Helps for maintaining correct joint position and muscle support during the function,

supports soft tissues around the elbow, provides proprioceptive effect and healing process by warming.


Features

Easy to fit and effective design, adjustable pressure by velcro system and removable pads.



4511 Manucare

Wrist support

cm	left/right
	XS : 14-15
	S : 15-16,5
	M : 16,5-18
	L : 18-19,5
	XL : 19,5-21
	XXL : 21-22,5

Indications

Tenosynovitis, carpal instabilities, CTS, after casting, arthralgia (reactive, rheumatic) gonglion pain, traumas of soft tissues of wrist .

Mode of Action

Stabilises the wrist, free range of motion for the thumb.


Features

High compliance and comfort thanks to D fabric, anatomically shaped and fitted design, reformable thermoplastic splint.



4531 Thenarcare

Thenar support

cm	left/right
	S : 10-14
	M : 14-18
	L : 18-22
	XL : 22-26

Indications

Degenerative changes of 1 st MCP, de quarvein, osteoarthritis, traumatic disorders

Mode of Action

Stabilises the MCP and CMP joints of thumb.


Features

High compliance and comfort thanks to D fabric, anatomically shaped and fitted design, easy to fit due to practical velcro system.



4521 Manucare plus

Wrist support with thumb

cm	left/right
	XS : 14-15
	S : 15-16,5
	M : 16,5-18
	L : 18-19,5
	XL : 19,5-21
	XXL : 21-22,5

Indications

MCP joint disorders, CMC joint degenerations, dequarvein, CTS + dequarvein.

Mode of Action

Stabilises the wrist and thumb.


Features

High compliance and comfort thanks to D fabric, anatomically shaped and fitted design, reformable thermoplastic splint.



4522 Manucare III plus

Wrist support with thumb

cm	left/right
	S : 12-14
	M : 14-16
	L : 16-18
	XL : 18-20
	XXL : 20-22

Indications

MCP joint disorders, CMC joint degenerations, dequarvein, CTS + dequarvein.

Mode of Action

Stabilises the wrist and thumb.


Features

High compliance and comfort due to anatomically shaped, easy to fit and remove sleeve design, practical velcro system and loops for pressure adjustment, aluminium splint.



4541 Manucare stable

Wrist support with stabilizing strap and wide thermoplastic splint

cm	left/right
	XS : 14-15
	S : 15-16,5
	M : 16,5-18
	L : 18-19,5
	XL : 19,5-21
	XXL : 21-22,5

Indications

Tenosynovitis, carpal instabilities, CTS, after casting, arthralgia (reactive, rheumatic) gonglion pain, traumas of soft tissues of wrist.

Mode of Action

Provides more stabilization with extra strap, free range of motion for the thumb.


Features

High compliance and comfort thanks to D fabric, anatomically shaped and fitted design, wide strap for pressure adjustment by velcro system, wide and reformable thermoplastic splint.



4542 Manucare III stable

Wrist support with stabilizing strap and wide thermoplastic splint

cm	left/right
	XS : 14-15
	S : 15-16,5
	M : 16,5-18
	L : 18-19,5
	XL : 19,5-21
	XXL : 21-22,5

Indications

Tenosynovitis, carpal instabilities, CTS, after casting, arthralgia (reactive, rheumatic) gonglion pain, traumas of soft tissues of wrist.

Mode of Action

Provides more stabilization with extra strap, free range of motion for the thumb.


Features

High compliance and comfort thanks to D fabric, anatomically shaped and fitted design, wide strap for pressure adjustment by velcro system, wide and reformable thermoplastic splint.



4911 Manucare easy

Knitted elastic wrist support

cm	
	
	S : 14-16
	M : 16-18
	L : 18-20
	XL : 20-22

Indications

Traumas of soft tissues of wrist, joint pain, tendinitis and arthritis, postraumatic and postoperative conditions, prophylactic

Mode of Action

Free range of motion for the

thumb, provides proprioceptive effect and healing process


Features

Anatomically knitted and supported with specially knitted zones at the wrist and hand, high elasticity, breathable comfort knitting system.



4320 Manucare active

Wrist support

cm	left/right
	
	XS : 13-14
	S : 14-16
	M : 16-18
	L : 18-20
	XL : 20-22

Indications

Overuse of wrist, after casting, tenosynovitis, carpal instabilities.

Mode of Action

Stabilises the wrist, free range of motion for the thumb.


Features

Compliance and comfort thanks to elastic material, anatomically shaped and fitted design, reformable thermoplastic splint, high quality of neoprene material.



4325 Manucare form

Wrist support with reformable thermoplastic splint

cm	left/right
	
	S : 14-16
	M : 16-18
	L : 18-20
	XL : 20-22

Indications

Tenosynovitis, carpal instabilities, CTS, after casting, arthralgia (reactive, rheumatic) ganglion pain, traumas of soft tissues of wrist

Mode of Action

Support both volar and dorsal of wrist, stabilises.

Features

Best compliance and comfort provided by reformable thermoplastic splint which is moulded to fit the particular patient's anatomy, easy to use due to practical velcro system, high quality of neoprene.



B10671 Strip splint

universal

B108 Frog splint

B109 Cot splint

B111 Baseball splint

S
M
L


B1070 Stax splints

1	5
2	5.5
3	6
4	7



4330 Thenarcare form

Thenar support with reformable thermoplastic splint

cm	left/right
	
	XS : 13-14
	S : 14-16
	M : 16-18
	L : 18-20
	XL : 20-22

Indications

Degenerative changes of 1st MCP, de quervein, osteoarthritis, traumatic disorders

Mode of Action

Stabilises the MCP and CMP joints of thumb.

Features


Best compliance and comfort provided by reformable thermoplastic splint which is moulded to fit the particular patient's anatomy, easy to use due to practical velcro system, high quality of neoprene material





5390 Thermal vest

Thermal vest

cm	
	L : 90-100 XL : 100-110

Indications


Decreased body stabilisation, insufficient blood circulation, to keep warm.

Mode of Action

Increases the blood flow, provides support and warmth.


Features

Anatomically shaped and fitted, high quality elastic and soft neoprene material.



3026 Shoulder support

Shoulder support

cm	
	S : 28-32 M : 32-35 L : 35-38

Indications

Humerus fractures, shoulder stabilisation.

Mode of Action

Stabilizes shoulder complex, prevents inferior subluxation/dislocation of glenohumeral joint.


Features

Comfortable and adjustable straps, easy to fit and remove due to practical velcro system, anatomic design covers the upper arm and shoulder.



5380 Thermal pant

Thermal pant

cm	
	S : 64-72 M : 72-80 L : 80-90 XL : 90-100 XXL : 100-110

Indications

Groin strains, hip abductors and flexor strains, stabilisation of the hip joint.

Mode of Action

Increases the blood flow, provides support and warmth.


Features

Anatomically shaped and fitted, high quality elastic and soft neoprene material.



5141 Groin support

Groin support

cm	left/right
	S : 49-51 M : 51-54 L : 54-57 XL : 57-60

Indications

Groin strains, hip abductors and flexor strains, stabilisation of the hip joint.

Mode of Action

Increases the blood flow and improves proprioception, lets easy flexing and stretching of muscles, provides support and warmth, relieves pressure.


Features

Anatomically shaped and fitted, upper width for better compliance, high quality neoprene material.



6040 Thigh support

Thigh support

cm	
	S : 48-50 M : 50-52 L : 52-54 XL : 54-56 XXL : 56-58

Indications

Hamstring and quadriceps strains, sprains and pulls, quadriceps contusions.

Mode of Action

Increases the blood flow and improves co-ordination and proprioception, provides support and warmth, relieves pressure.


Features

Anatomically shaped and fitted, upper width for better compliance, high quality neoprene material.



6060 Calf support

Calf support

cm	
	S : 31-33 M : 33-35 L : 35-37 XL : 37-39 XXL : 39-41

Indications

Compartment syndrome, contusion, shin splints, sprain, muscle strains.

Mode of Action

Increases the blood flow and improves co-ordination and proprioception, provides support and warmth, relieves pressure.

Features

Anatomically shaped and fitted, upper width for better compliance, high quality neoprene material.



6157 Genucare jumper's knee

Knee support for jumpers knee

cm	
S	: 30-35
M	: 35-40
L	: 40-45

Indications

Osgood schlatter jumper's knee.

Mode of Action

Relief the pain due to applied pressure on patellar tendon, stabilises patella due straps.

Features

Anatomically shaped and fitted, supported with straps, high quality neoprene material.



6158 Genucare Luxa

Knee support with lateral pad

cm	
S	: 30-35
M	: 35-40
L	: 40-45

Indications

Luxation of the patella, patello-femoral malalignments, after lateral release surgery.

Mode of Action

Stabilises the knee cap laterally due to pad and two straps, unloads the patella tendon.

Features

Anatomically shaped and fitted, supported with patella pad and straps, high quality neoprene.



6875 Knee strap Inno

Patellar tendon strap

cm	
1	: 25-35
2	: 35-45

Indications

Pain due to patello-femoral malalignments, mild instabilities of knee joint.

Mode of Action

Supports the knee joint , relieves proprioception.

Features

Anatomically shaped fitting design which is ideal for sports and exercise, compressive foam pad and strap allowing to adjust the compression, durable loop fabric.



6876 Patellar tendon support

Patellar tendon strap

universal

Indications

Pain due to patello-femoral malalignments, mild instabilities of knee joint.

Mode of Action

Supports the knee joint , relieves proprioception.

Features

Anatomically shaped, supported with patella pad, high quality neopren material.



6130 Genucare ligament

Knee support with patella pad and soft splints

cm	
XS	: 27-30
S	: 30-33
M	: 33-36
L	: 36-39
XL	: 39-42
XXL	: 42-45

Indications

Arthritis, osteoarthritis, irritation posttraumatic and post-operative, tendomyopathies, feeling of instability.

Mode of Action

Supports and stabilises the kneecap and joint, patella ring and patella cut out give extra support, improves proprioception.

Features

High compliance due to wrap around velcro straps, medial and lateral spiral stays, high quality neoprene material.



6135 Genucare ligament wrap

Knee support with soft splints

cm	
S	: 30-33
M	: 33-36
L	: 36-39
XL	: 39-42
XXL	: 42-45

Indications

Arthritis, osteoarthritis, irritation posttraumatic and post-operative, tendomyopathies, feeling of instability.

Mode of Action

Knee support gives extra support and stability to the kneecap and joint, patella ring and patella cut out give extra support, improves proprioception.

Features

High compliance due to wrap around version and practical velcro closing system at the front, comfortable wearing due to special popliteal zone material, medial and lateral spiral stays, high quality neoprene material.



6140 Genucare basic open

Knee support open patella

cm	
	S : 30-34 M : 34-38 L : 38-42 XL : 42-46 XXL : 46-50

Indications

Osteoarthritis, efusion, other painful conditions of knee.

Mode of Action

Increases the blood flow and improves co-ordination and proprioception, provides support and warmth, relieves pressure.

Features

Anatomically shaped and fitted, high quality neoprene material.



6141 Genucare patella

Knee support with patella pad

cm	
	S : 30-33 M : 33-36 L : 36-39 XL : 39-42 XXL : 42-45

Indications

Osteoarthritis, efusion, other painful conditions of knee.

Mode of Action

Knee support with pad that gives extra support and stability to the kneecap and joint, warmth, improves proprioception.

Features

Anatomically shaped and fitted, supported with patella ring, high quality neoprene material.



6150 Genucare stable

Knee support with patella pad and hinge

cm	
	S : 30-34 M : 34-38 L : 38-42 XL : 42-46 XXL : 46-50

Indications

ACL and PCL injuries, 1st and 2nd grade meniscus degenerations, collateral ligament injury, arthritis (progressive rheumatoid arthritis).

Mode of Action

Supports and stabilises the knee and joint, patella ring and patella cut out give extra support, improves

proprioception, prevents hyperextension.

Features

High compliance due to wrap around velcro straps, medial and lateral polycentric hinges, high quality neoprene material.



6155 Genucare stable wrap

Knee support with hinge

cm	
	S : 30-34 M : 34-38 L : 38-42 XL : 42-46 XXL : 46-50

Indications

ACL and PCL injuries, 1 st and 2 nd grade meniscus degenerations, collateral ligament injury, arthritis (progressive rheumatoid arthritis).

Mode of Action

Supports and stabilises the knee and joint, patella ring and patella cut out give extra support, improves

proprioception, prevents hyperextension.

Features

High compliance due to wrap around version and practical velcro closing system at the front, comfortable wearing due to special popliteal zone material, medial and lateral polycentric hinges, high quality neoprene material.



6120 Genucare basic

Knee support

cm	
	S : 30-34 M : 34-38 L : 38-42 XL : 42-46 XXL : 46-50

Indications

Osteoarthritis, efusion, other painful conditions of knee.

Mode of Action

Increases the blood flow and improves co-ordination and proprioception, provides support and warmth, relieves pressure.

Features

Anatomically shaped and fitted, high quality neoprene material.



6160 Genucare hyperX

Knee brace for prevention of hyper extension

cm	
	S : 30-33 M : 33-36 L : 36-39 XL : 39-42 XXL : 42-45

Indications

ACL and PCL injuries, 1st and 2nd grade meniscus degenerations, collateral ligament injury, arthritis (progressive rheumatoid arthritis).

Mode of Action

Supports and stabilises the knee and joint, patella ring and patella cut out give extra support, improves proprioception, prevents hyperextension.

Features

Prevention of hyperextension thanks to crossed straps at the back, high compliance due to velcro straps at the front, easy to pull on and comfortable wearing due to special popliteal zone material, medial and lateral polycentric hinges, high quality neoprene material.



6510 Knee support wool

Knee support

cm	
S	: 30-33
M	: 33-36
L	: 36-39
XL	: 39-42
XXL	: 42-45

Indications

Rheumatism, knee pain, to keep knee warm.

Mode of Action

Increases the blood flow and improves co-ordination provides support and warmth, relieves proprioception.

Features

High quality and elastic angoré.



6911 Genucare Easy

Knitted elastic knee support

	A	B
S	: 43-46	31-34
M	: 46-49	34-37
L	: 49-52	37-40
XL	: 52-55	40-43
XXL	: 55-58	43-46

Indications

All inflammations of knee, posttraumatic and post-operative, conditions, arthritis, osteoarthritis, tendomyopathies, feeling of instability.

Mode of Action

Supports the kneecap and joint, knitted elastic material

Features

Breathable knitting system, high compliance with specially knitted zones at the popliteal area, high elasticity.



7620 Malleocare basic

Ankle support

cm	
S	: 20-22
M	: 22-24
L	: 24-26
XL	: 26-28
XXL	: 28-30

Indications

Tendinitis, joint effusions, arthritis and osteoarthritis, ligament weaknesses, for compression.

Mode of Action

Increases the blood flow and reduces stiffness. Flexible

support that fits easily in shoes.

Features

Anatomically shaped and fitted, easy to fit and remove due to practical velcro system, high quality neoprene material.



7622 Malleocare basic plus

Ankle support with strap

cm	
S	: 20-22
M	: 22-24
L	: 24-26
XL	: 26-28
XXL	: 28-30

Indications

Tendinitis, joint effusions, arthritis and osteoarthritis, ligament weaknesses, for compression.

Mode of Action

Increases the blood flow and reduces stiffness, stabilizes the ankle.

Features

Anatomically shaped and fitted, easy to fit and remove due to practical velcro system, high quality neoprene material.



7530 Ankle support light

Ankle support

cm	
XS	: 17-20
S	: 20-22
M	: 22-24
L	: 24-26
XL	: 26-28
XXL	: 28-30

Indications

Post operative and posttraumatic irritation, to support joint, painful conditions of ankle.

Mode of Action

Elastic material helps circulation and increases the

blood flow, reduces stiffness, improves proprioception.

Features

Light and flexible, easily fits in shoes.



7630 Malleocare flex

Ankle support with pads

cm	
XS	: 17-20
S	: 20-22
M	: 22-24
L	: 24-26
XL	: 26-28
XXL	: 28-30

Indications

Tendinitis, joint effusions, arthritis and osteoarthritis, ligament weaknesses, for compression.

Mode of Action

The pad has a massage effect which helps circulation and

increases the blood flow, reduces stiffness and stabilises the ankle.

Features

Anatomically shaped and fitted, easy to fit and remove due to practical velcro system, high quality neoprene material.



7640 Malleocare stable

Ankle brace with plastic splint

cm	
	S : 20-22 M : 22-24 L : 24-26 XL : 26-28 XXL : 28-30

Indications

Stable fractures of malleolus, chronic instability of ankle, ligament injuries, post-operative care, joint effusion/swelling.

Mode of Action

Protects the ankle, let to move and stabilises.

Features

Plastic insert guarding malleolus, anatomically shaped and fitted, easy to fit and remove due to practical velcro system, high quality neoprene material.



7642 Malleocare stable plus

Ankle brace with stabilising strap

cm	
	S : 20-22 M : 22-24 L : 24-26 XL : 26-28 XXL : 28-30

Indications

Stable fractures of malleolus, chronic instability of ankle, ligament injuries, post-operative care, joint effusion/swelling.

Mode of Action

Protects the ankle, let to move and stabilises.

Features

Plastic insert guarding malleolus, anatomically shaped and fitted, easy to fit and remove due to practical velcro system, 8 bandage for mor stabilisation high quality neoprene material.



7675 Malleocare form

Ankle brace with reformable thermoplastic splint

cm	
	S : 20-22 M : 22-24 L : 24-26 XL : 26-28 XXL : 28-30

Indications

ACL and PCL injuries, 1st and 2nd grade meniscus degenerations, collateral ligament injury, arthritis (progressive rheumatoid arthritis).

Mode of Action

Supports and stabilises the knee and joint, patella ring and patella cut out give extra support, improves proprioception, prevents hyperextension.

Features

Prevention of hyperextension thanks to crossed straps at the back, high compliance due to velcro straps at the front, easy to pull on and comfortable wearing due to special popliteal zone material, medial and lateral polycentric hinges, high quality neoprene material.



7911 Malleocare Easy

Knitted elastic ankle support

cm	
	S : 20-22 M : 22-24 L : 24-26 XL : 26-28 XXL : 28-30

Indications

Chronic and acute inflammation of ankle, ligament injuries, prophylactic purposes for sport.

Mode of Action

The elastic comfort knitting provides compression which helps circulation and increases the blood flow, reduces

stiffness and stabilise, improves proprioception.

Features

Breathable comfort knitting system, anatomically knitted, high compliance with specially knitted zones at the ankle bending area.



7110 Silicon insole

shoe size	
S	: 36-38
M	: 39-40
L	: 41-42
XL	: 43-44
XXL	: 45-46

Indications

Pes planus, heel pain, pressure redistribution, plantar fasciitis, fat pad, arthralgia.

Mode of Action

Helps to solve the cause of foot pain.

Features

Shock absorbing blue areas, air ventilation system, anatomically shaped structure with arch support, perfect silicon material.



7120 Heel cup

shoe size	
S	: 36-38
M	: 39-41
L	: 42-44

Features

Shock absorbing blue areas, anatomically shaped structure, perfect silicon material.

Indications

Achillodynia- tendomyopathy, calcaneal spur, haglund's deformity, plantar fasciitis, calcaneal pain.

Mode of Action

Provide a soft tissue around the calcaneus and relieve pressure on heel spur.



7130 Metatarsal cushion

universal

Indications

Conditions to support transverse arch, metatarsalgia, irritations of metatarsal area, collozites of metatars head.

Mode of Action

Silicon cushion provide a comfortable supports for the transvers arch.

Features

Comfortable and easy use non-slip knitted material, anatomically shaped and placed silicon support.



7140 Toe separator

shoe size	
S	: 36-38
M	: 39-41
L	: 42-44

Features

Comfortable and easy use non-slip silicon support.

Indications

Hallux valgus - conservative and post-operative treatment.

Mode of Action

Unload the joints of big toe with a tension and helps moving the big toe to the right position



7310 Insole 3/4 stable

Insole with steel support

shoe size	
0	: 33-34
1	: 35-36
2	: 37-38
3	: 39-40
4	: 41-42
5	: 43-44
6	: 45-46

Indications

Pes planus, pes cavus, metatarsalgia.

Mode of Action

Supports and maintain normal longitudinal arch and provides cushioning and support for all metatarsal heads in normal transverse arch position.

Features

Anatomically shaped, made of unique leather material and supported with stainless steel. A soft cushion is placed in the metatars area for more comfort. Durable and long life material.



7320 Insole 3/4

Insole

shoe size	
35	40
36	41
37	42
38	43
39	44

Features

Anatomically shaped, made of artificial leather material. A soft cushion is placed in the metatarsal area for comfort.

Indications

Pes planus, pes cavus, metatarsalgia.

Mode of Action

Supports and maintain normal longitudinal arch and provides cushioning and support for all metatarsal heads in normal transverse arch position.



1912 Back pillow
Lumbar back pillow

universal

Features

Anatomically shaped.
The knitted cover is made of polypropylene. It breathes and provides moisture transfer.



1913 Sit pillow
Sit pillow

universal

Features

Anatomically shaped.
The knitted cover is made of polypropylene. It breathes and provides moisture transfer.



1911 Memory foam pillow
Memory foam neck pillow

universal

Features

Anatomically shaped, comfortable and soft. Made of memory foam visco.
The knitted cover breathes, provides moisture transfer.

SMOOTH PILLOWS



orthocare®

New products with
authentic designs

*ROM elbow and knee orthoses
Ankle foot orthoses*



Arden Medikal factory, Ankara



Please scan
the code
for complete
product
catalogue.

©    in
orthocaretr

ARDEN MEDİKAL

Saray Mahallesi 100. Yıl Bulvarı
No:16 K.Kazan Ankara
06980 Türkiye
+90 312 394 14 54
www.ardenmedikal.com.tr
export@ardenmedikal.com.tr